



北京大学
PEKING UNIVERSITY



北京核磁共振中心
Beijing NMR Center, Peking University

核磁共振上机培训

李红卫

2020-10



北京大学
PEKING UNIVERSITY



北京核磁共振中心
Beijing NMR Center, Peking University

核磁共振上机基础

李红卫

2020-10-25

基本原理


基本条件：
稳定磁场 B_0



核 磁 共 振

检测对象：自旋量子数 I 不为0

自旋量子数(I)	示例
半整数 (1/2, 3/2...)	^1H , ^{13}C , ^{15}N
整数 (1, 2, ...)	^2D , ^{14}N



检测内容：
原子核进动频率


$$\nu = \frac{\gamma B_0}{2\pi}$$

γ : 旋磁比, 与不同原子核的性质相关;

ν : 进动频率; B_0 : 磁场强度;

➤ Spin, Precession and Larmor frequency

化学位移及分辨率




$$\nu_{real} = \frac{\gamma B_{eff}}{2\pi}$$


$$\delta(ppm) = \frac{\nu - \nu_{ref}}{\nu_{BF1}} \times 10^6$$

场强越高，分辨率越高

核磁共振信号产生




Bulk Magnetization




1D experiments: recycle delay pulses detection



核磁共振仪器及注意事项




磁体




Attention

- 使用心脏起搏器和体内有金属植入者禁止进入核磁共振实验室。
- 严禁携带以下物品进入核磁共振仪器实验室：铁磁性物质（铁制工具，钥匙等），电子设备（手机，手表等），电子卡片（银行卡等），除核磁样品以外的化学品，食品和饮料。
- 磁体处于悬浮状态，上样时禁止触碰磁体

2K磁体





探头

进入核磁实验室对仪器的任何操作都需要具有核磁上机权限的人陪同

Attention

- 使用完好无损的核磁管
- 无自动上样器的仪器，上样前需确认Lift气流
- 上样前要将核磁管擦拭干净，特别是刚从低温取出的样品
- 使用合适的脉冲功率







控制器



信号通路



实验操作流程



Sample preparation

Field preparation (lock, atmm, topshim)

Pulse calibration

Experiment-related parameter



样品准备

- 样品的制备

- 非标记样品； ^{15}N 标记样品； ^{15}N 和 ^{13}C 标记样品；其他特殊标记样品；
- 样品浓度：~mM量级；
- 溶剂的体积为450-550 μl ；
- 选择合适的缓冲体系， $\text{pH} < 7.5$ ；DSS； D_2O ；
- 核磁管一般为5mm内径的核磁管；

- 进样与出样

- Bruker 自动进样系统
- Lift命令进样







磁场准备

Δ Solvent	Description
Acetic	acetic acid-d4
Acetone	acetone-d6
C6D6	benzene-d6
CD2Cl2	dichloromethane-d2
CD3CN	acetonitrile-d3
CD3CN_SPE	LC-SPE Solvent (Acetonitrile)
CD3OD_SPE	LC-SPE Solvent (Methanol-d4)
CDCl3	chloroform-d
CH3CN+D2O	HPLC Solvent (Acetonitrile/D2O)
CH3OH+D2O	HPLC Solvent (Methanol/D2O)
D2O	deuteriumoxide
D2O_salt	deuteriumoxide with salt
Dioxane	dioxane-d8
DMF	N,N-dimethylformamide-d7
DMSO	dimethylsulfoxide-d6
EtOD	ethanol-d6
H2O+D2O	90%H2O and 10%D2O
H2O+D2O_salt	90%H2O and 10%D2O with salt
HDMSO	90%DMSO and 10%DMSO-d6
Juice	fruit juice
MeOD	methanol-d4
Plasma	blood plasma
Pyr	pyridine-d6
T_H2O+D2O+Me4NCl	(CD3)4NCl in 90%H2O and 10%D2O, for NMR thermometer
T_H2O+D2O+NaAc	sodium acetate in 90%H2O and 10%D2O, for NMR thermometer
T_H2O+D2O+Pivalate	pivalate-d9 in 90% H2O and 10% D2O, for NMR thermometer
T_MeOD	methanol-d4, for NMR thermometer
TFE	trifluoroethanol-d3
THF	tetrahydrofuran-d8
Tol	toluene-d8
Urine	urine



lock



topshim

^1H 的刻度

利用 360° 信号最弱的原理来刻度 90° 脉冲长度: p1和plw1





^1H 的刻度实验

```
;zg
;avance-version (12/01/11)
;1D sequence
;
;$CLASS=HighRes
;$DIM=1D
;$TYPE=
;$SUBTYPE=
;$COMMENT=
```


```
1 ze
2 30m
d1 recycle delay
p1 ph1 pulses
go=2 ph31
30m mc #0 to 2 F0(zd) detection
exit
```

```
#include <Avance.incl>
```

```
ph1=0 2 2 0 1 3 3 1
ph31=0 2 2 0 1 3 3 1
```


```
"acqt0=-p1*2/3.1416"
```

```
;p1 : f1 channel - power level for pulse (default)
;p1 : f1 channel - high power pulse
;d1 : relaxation delay; 1-5 * T1
;ns: 1 * n, total number of scans: NS * TD0
;$ld: zg,v 1.10.8.1 2012/01/31 17:56:41 ber
Exp $
```



^1H 的刻度实验具体参数

Spectrum		ProcPars		AcquPars		Title		PulseProg		Peaks		Integrals		Sample		Structure		Plot		Fid	
Probe: not defined																					
General Channel f1																					
General																					
PULPROG	zg			...		E		Pulse program for acquisition													
TD	65536			Time domain size																	
SWH [Hz, ppm]	8012.82			20.0254			Sweep width														
AQ [sec]	4.0894465			Acquisition time																	
RG	16			Receiver gain																	
DW [μsec]	62.400			Dwell time																	
DE [μsec]	10.00			Pre-scan-delay																	
D1 [sec]	1.00000000			Relaxation delay; 1-5 * T1																	
DS	0			Number of dummy scans																	
NS	1			1 * n, total number of scans: NS * TD0																	
TD0	1			Dimension of accumulation loop																	
Channel f1																					
SFO1 [MHz]	400.1318700			Frequency of ch. 1																	
O1 [Hz, ppm]	1869.97			4.673			Frequency of ch. 1														
NUC1	1H			Edi...			Nucleus for channel 1														
P1 [μsec]	8.61			F1 channel - high power pulse																	
PLW1 [W, -dBW]	12.589			-11.00			F1 channel - power level for pulse (default)														





北京大学
PEKING UNIVERSITY



北京核磁共振中心
Beijing NMR Center, Peking University

上机练习核磁实验基本流程和 ^1H 脉冲刻度



北京大学
PEKING UNIVERSITY



北京核磁共振中心
Beijing NMR Center, Peking University


水峰压制的 ^1H 谱

李红卫


2020-11-01




水峰压制的意义



不同脉冲及激发范围差异



脉冲功率和时间的不同导致其激发的化学位移范围不同





edprosol

File Edit View Help

Modified Observe and Saved Decouple Prosol Parameter Set for:

Probe: Z118687_0006 CP TCI 500S2 H-C/N-D-05 Z

Select ...

Solvent: generic

Observe

Decouple

1H

Nucleus

13C

Observe

Decouple

Observe Comment: Default 1H obs 500

Decouple Comment: Default 13C dec 500

90 deg. Pulses

HR Square Pulses

HR Shape Pulses

Others

Observe

Decouple

Nucleus	Pulse Width[μs]	Att. Lvl.[dB]	Set	Pulse Width[μs]	Att. Lvl.[dB]	Set	Nucleus
1H	11.45	-8.26		8.00	-8.26		1H
2H	67.80	-16.23		67.80	-16.23		2H
13C	12.00	-19.68		12.00	-19.68		13C
15N	25.00	-21.52		25.00	-21.52		15N

Observe

Decouple

Nucleus	Pulse Width[μs]	Att. Lvl.[dB]	Set	Pulse Width[μs]	Att. Lvl.[dB]	Set	Nucleus
1H	11.45	-8.26		8.00	-8.26		1H
2H	67.80	-16.23		67.80	-16.23		2H
13C	12.00	-19.68		12.00	-19.68		13C
15N	25.00	-21.52		25.00	-21.52		15N

Observe

Decouple

Nucleus	Pulse Width[μs]	Att. Lvl.[dB]	Set	Pulse Width[μs]	Att. Lvl.[dB]	Set	Nucleus
1H	11.45	-8.26		8.00	-8.26		1H
2H	67.80	-16.23		67.80	-16.23		2H
13C	12.00	-19.68		12.00	-19.68		13C
15N	25.00	-21.52		25.00	-21.52		15N

Observe

Decouple

Nucleus	Pulse Width[μs]	Att. Lvl.[dB]	Set	Pulse Width[μs]	Att. Lvl.[dB]	Set	Nucleus
1H	11.45	-8.26		8.00	-8.26		1H
2H	67.80	-16.23		67.80	-16.23		2H
13C	12.00	-19.68		12.00	-19.68		13C
15N	25.00	-21.52		25.00	-21.52		15N

Last Save

Print

Copy to Solvent

Copy to Probe

Save

edprosol

22




Observe			Decouple		
PuW[μs]	A[dB]	#	PuW[μs]	A[dB]	
cpd	45.00	3.63	0	cpd	75.00 -3.71
TOCSY spinlock	25.00	-1.48	1	TOCSY spinlock	25.00 -13.31
ROESY spinlock (cw, RF field)	110.00	11.39	2	ROESY spinlock	125.00 0.67
presat. (cw irradiation, RF field)	5000.00	44.54	3	presat. (cw irradiation, RF field)	250.00 6.69
			4	2nd cpd (power gated)	98.70 -1.38
			5	low power cpd	140.00 1.65
			6	bilev cpd (cw part)	35.00 -10.39
TOCSY/hetero T. (med. selectivity)	52.00	4.88	7	TOCSY/hetero T. (med. selectivity)	52.00 -6.95
TOCSY/hetero T. (high selectivity)	78.00	8.40	8	TOCSY/hetero T. (high selectivity)	78.00 -3.43
			9	TOCSY/hetero T. (very high selectivity)	250.00 6.69
cleanex spinlock	52.00	4.88	10		
ROESY pulsed (90°)	108.00	11.23	11		
low power presat. (cw irradiation, RF field)	25000.00	58.52	12	low power presat. (cw irradiation, RF field)	
			13		
			14	homodecoupling	750.00 16.23
			15		
			16		
			17	low power cpd (2)	480.00 12.36
			18	180° (matched field for Calp)	128.90 -5.08
			19		
			20		
			21		
			22		
			23		

Observe				Decouple			
Filename	PuW[μs]	A[dB]	#	Filename	PuW[μs]	A[dB]	
selective excitation	Gaus1_270.1000	80000.00	51.37	0	selective excitation	Q5.1000	40000.00 25.50
selective refocussing	Gaus1_180r.1000	80000.00	54.89	1	select. inversion/refocussing	Q3_surbop.1	40000.00 26.37
bandsel. excitation	Q5.1000	10000.00	25.29	2	bandsel. excitation	Q5.1000	3600.00 4.59
bandsel. inv./refoc.	Q3.1000	10000.00	28.15	3	bandsel. inv./refoc.	Q3_surbop.1	3792.00 5.91
off-resonance presat. (power)	Squa100.1000	100000.00	70.56	4	adiabatic inversion	Crp60,0.5,20.1	500.00 -13.11
90° flip back (H2O)	Squa100.1000	1000.00	30.56	5	adiabatic refocussing	Crp60comp.4	2000.00 -13.11
2nd 90° flip back (H2O)	Sinc1.1000	4000.00	38.01	6	Cal/CO 90°	Q5.1000	384.00 -14.85
90° WET	Sinc1.1000	20000.00	51.98	7	Cal/CO 90° timerev.	Q5tr.1000	384.00 -14.85
120° NH region	Pc9_4_120.1000	2880.00	19.19	8	Cal/CO 180°	Q3.1000	307.00 -13.94
180° NH region I	Rsnob.1000	960.00	10.78	9	Calpha sel. 90°	Q5.1000	1800.00 -1.43
90° NH region I	Pc9_4_90.1000	2640.00	20.93	10	Calpha sel. 90° timerev.	Q5tr.1000	1800.00 -1.43
90° NH region I timerev.	Pc9_4_90.1000	2640.00	20.93	11	Calpha sel. 180°	Q3.1000	1000.00 -3.68
180° NH region II	Reburp.1000	1680.00	7.09	12	adiabatic decoupling	Crp32,1.5,20.2	1500.00 -3.39
90° NH region II	Eburp2.1000	2040.00	12.47	13	adiab. decoupling (bilev par)	Crp32,1.5,20.2	750.00 -9.41
90° NH region II timerev.	Eburp2tr.1000	2040.00	12.47	14	180° short broadband	Bip720,100,10.1	192.00 -13.66
180° broadband	Bip720,50,20.1	240.00	0.11	15	180° medium selectivity	Q3_surbop.1	888.00 -6.70
cleanex 90° (H2O)	Gaus1_90.1000	5000.00	36.83	16	180° high selectivity	Q3_surbop.1	2520.00 2.36
cleanex 180° (H2O)	Gaus1_180r.1000	7500.00	34.33	17	90° high selectivity	Q5.1000	4440.00 6.41
z-spoil (adiabatic)	Crp60,20,20.10	20000.00	14.74	18	90° high selectivity timerev.	Q5tr.1000	4440.00 6.41
180° H1' (na_)	Reburp.1000	3960.00	14.54	19	x-filter (adiabatic 180)	Crp60_xfilt.2	1895.00 -6.93
				20	Ca or CO decoupling	Q3_surbop.1	972.00 -5.92
180° H2O (selective)	Gaus1_180i.1000	120000.00	58.42	21	simult. Ca + CO selective	Q3Ca_CaCO.1000	922.00 -10.41
180° Halpha sel.	Reburp.1000	1920.00	8.25	22	Cbeta decoupling		
90° Hmethyl sel.	Pc9_4_90.1000	7200.00	29.65	23	na_ 90°	Q5.1000	1200.00 -4.96



presat压水 (zgpr)






PULPROG	zgpr	...	E	Pulse program for acquisition
TD	32768			Time domain size
SWH [Hz, ppm]	8012.82	16.0214		Sweep width
AQ [sec]	2.0447233			Acquisition time
RG	1			Receiver gain
DW [μsec]	62.400			Dwell time
DE [μsec]	10.00			Pre-scan-delay
D1 [sec]	1.000000000			Relaxation delay; 1-5 * T1
d12 [sec]	0.00002000			Delay for power switching [20 usec]
DS	0			Number of dummy scans
NS	1			1 * n, total number of scans: NS * TD0
TD0	1			Number of averages in 1D
Channel f1				
SFO1 [MHz]	500.1323366			Frequency of ch. 1
O1 [Hz, ppm]	2336.61	4.672		Frequency of ch. 1
NUC1	1H	Edit...		Nucleus for channel 1
P1 [μsec]	8.000			F1 channel - 90 degree high power pulse
PLW1 [W, dB]	6.7	8.26		F1 channel - power level for pulse (default)
PLW9 [W, dB]	1.7152e-005	47.66		F1 channel - power level for presaturation

预饱和压水时间

预饱和压水功率

脉冲梯度场




$$\nu = \frac{\gamma B_0}{2\pi}$$

消除不需要的信号



watergate压水 (zgpgpwg)





PULPROG	zgpgwg	E	Pulse program for acquisition	
TD	32768	Channel f1		
SWH [Hz, ppm]	8012.82	16.0214		
AQ [sec]	2.0447233	SFO1 [MHz] 500.1323506 Frequency of ch. 1		
RG	203	O1 [Hz, ppm] 2350.61 4.700 Frequency of ch. 1		
DW [μsec]	62.400	NUC1 1H Edit... Nucleus for channel 1		
DE [μsec]	10.00	P1 [μsec] 8.000 F1 channel - 90 degree high power pulse		
D1 [sec]	2.000000000	p2 [μsec] 16.00 F1 channel - 180 degree high power pulse		
d12 [sec]	0.00002000	P11 [μsec] 1000.000 F1 channel - 90 degree shaped pulse		
D16 [sec]	0.000200000	PLW0 [W, dB] 0 1000.00 0W		
DS	4	PLW1 [W, dB] 6.7 -8.26 F1 channel - power level for pulse (default)		
NS	8	SPNAM 1 Sinc1.1000 ... E File name for SP1		
TD0	1	SPOAL1 0.500 Phase alignment of freq. offset in SP1		
			SPOFFS1 [Hz] 0 Offset frequency for SP1	
			SPW1 [W, -dBW] 0.0004288 33.68 F1 channel - shaped pulse 90 degree	
Gradient channel				
GPNAM 1 SMSQ10.100 ... E SMSQ10.100			20% → 压水梯度场	
GPZ1 [%] 20.00			Homospoil/gradient pulse	
P16 [μsec] 1000.000				


watergate

水的选择性激发

压水梯度场



组合脉冲压水watergate w5 pulse





PULPROG	zgpgw5	...	E
TD	32768		
SWH [Hz, ppm]	8012.82	16.0214	
AQ [sec]	2.0447233		
RG	203		
DW [μsec]	62.400		
DE [μsec]	10.00		
D1 [sec]	2.000000000		
D16 [sec]	0.000200000		
D19 [sec]	0		
DS	4		
NS	8		
TD0	1		

组合脉冲时间间隔

Pulse program for acquisition
Time domain size
Sweep width
Acquisition time
Receiver gain
Dwell time
Pre-scan-delay

Channel f1

SFO1 [MHz]	500.1323506	
O1 [Hz, ppm]	2350.61	4.700
NUC1	1H	Edit...
P1 [μsec]	8.000	
P27 [μsec]	8.000	
PLW1 [W, dB]	6.7	-8.26
PLW18 [W, dB]	6.7	-8.26

组合脉冲脉宽和功率

Frequency of ch. 1
Frequency of ch. 1
Nucleus for channel 1
F1 channel - 90 degree high power pulse
F1 channel - 90 degree pulse at p18
F1 channel - power level for pulse (default)
F1 channel - power level for 3-9-19-pulse (watergate)

Gradient channel

GPnam 1	SMSQ10.100	...	E	SMSQ10.100
GPZ1 [%]	34.00			34%
GPnam 2	SMSQ10.100	...	E	SMSQ10.100
GPZ2 [%]	22.00			22%
P16 [μsec]	1000.000			Homospoil/gradient pulse

压水梯度场

watergate W5



北京大学
PEKING UNIVERSITY



北京核磁共振中心
Beijing NMR Center, Peking University

上机练习压水 ^1H 谱 zgpgpw5和zgpr



2D谱.....

- J coupling
- T1 and T2
- INEPT
- NQSQC



北京大学
PEKING UNIVERSITY



北京核磁共振中心
Beijing NMR Center, Peking University

上机练习 ^1H - ^{15}N HSQC



北京大学
PEKING UNIVERSITY



北京核磁共振中心
Beijing NMR Center, Peking University

3D三共振实验

李红卫

2020-11-15






脉冲构成基本模块


- Hard pulse
- Shape pulse
- Composed pulse
- Gradient
- Delay and coupling
- Phase cycle





1D-2D-3D实验



+




whatever nucleus X (ppm)









1D experiments:




2D experiments:




1D: ~min;
2D: ~hours;
3D: ~days



⋮





3D experiments:



Multidimension: Cycles= $n1*n2*.....$

J-coupling:



- 1、磁矩传递
- 2、化学位移标记
- 3、化学位移标记时去耦
- 4、去除水和其他不需要的信号



General

PULPROG	hncagpwg3d	d10 [sec]	0.00584650
TD	1024	d11 [sec]	0.03000000
SWH [Hz, ppm]	6009.62 12.0161	D16 [sec]	0.000200000
AQ [sec]	0.0851968	d21 [sec]	0.00550000
RG	512	d23 [sec]	0.01200000
DW [μsec]	83.200	d26 [sec]	0.00230000
DE [μsec]	20.00	d29 [sec]	0.00029750
d0 [sec]	0.00000300	d30 [sec]	0.00584650
D1 [sec]	1.100000024	DELTA [sec]	0.00000600
		DELTA1 [sec]	0.00645500
		DELTA2 [sec]	0.00008800
		DS	32
		in0 [sec]	0.00012420
		in10 [sec]	0.00018970
		in29 [sec]	0.00018970
		in30 [sec]	0.00018970
		INF1 [μsec]	248.40
		INF2 [μsec]	758.80
		NS	32
		TDav	0

incremented delay (F2 in 3D) = $d_{23}/2 - p_{14}/2$
 Delay for disk I/O [30 msec]
 Delay for homospoil/gradient recovery
 $1/(2J(NH))$ [5.5 msec]
 $1/(4J(NCa))$ [12 msec]
 $1/(4J'(NH))$ [2.3 msec]
 incremented delay (F2 in 3D) = $d_{23}/2 - p_{14}/2 - p_{26} - d_{21} - 4u$
 decremented delay (F2 in 3D) = $d_{23}/2 - p_{14}/2$
 $DELTA = d_0^2 + \text{larger}(p_{14}, p_{22}) - p_{14}$
 $DELTA1 = d_{23} - d_{21} - p_{26}$
 $DELTA2 = d_{26} - p_{16} - p_{11} - 12u$
 ≥ 16
 $1/(2 * SW(Ca)) = DW(Ca)$
 $1/(4 * SW(N)) = (1/2) DW(N)$
 $= in_{10}$
 $= in_{10}$
 $1/SW(Ca) = 2 * DW(Ca)$
 $1/SW(N) = 2 * DW(N)$
 $8 * n$
 Number of averages in nD

General Parameter



Channel f1

SFO1 [MHz]	500.1323506		Frequency of ch. 1
O1 [Hz, ppm]	2350.61	4.700	Frequency of ch. 1
NUC1	1H	Edit...	Nucleus for channel 1
CPDPRG 1	waltz65	... E	File name for cpd1
P1 [μsec]	11.700		F1 channel - 90 degree high power pulse
p2 [μsec]	23.40		F1 channel - 180 degree high power pulse
P11 [μsec]	1000.000		F1 channel - 90 degree shaped pulse [1 msec]
P26 [μsec]	45.000		F1 channel - 90 degree pulse at pl19
PCPD1 [μsec]	45.00		F1 channel - 90 degree pulse for decoupling sequence
PLW0 [W, dB]	0	1000.00	0W
PLW1 [W, dB]	6.2	-7.92	F1 channel - power level for pulse (default)
PLW19 [W, dB]	0.4	3.98	F1 channel - power level for CPD/BB decoupling
SPNAM 1	Squa100.1000	... E	File name for SP1
SPOAL1	0.500		Phase alignment of freq. offset in SP1
SPOFFS1 [Hz]	0		Offset frequency for SP1
SPW1 [W, -dBW]	0.00079433	31.00	F1 channel - shaped pulse 90 degree (H2O on resonance)

1H Parameter:
45us ~ plw9



Channel f2

SFO2 [MHz]	125.7645794	Frequency of ch. 2	
O2 [Hz, ppm]	6790.92	54.000	Frequency of ch. 2
NUC2	13C	Edit...	Nucleus for channel 2
CNST21	174.0000000	CO chemical shift (offset, in ppm)	
P13 [μsec]	384.000	F2 channel - 90 degree shaped pulse	
P14 [μsec]	307.000	F2 channel - 180 degree shaped pulse	
PLW2 [W, dB]	0	1000.00	Power PLW2
SPNAM 2	Q5.1000	E	File name for SP2
SPOAL2	1.000	Phase alignment of freq. offset in SP2	
spoffs2 [Hz]	0	spoffs2=0	
SPW2 [W, -dBW]	26.95	-14.31	F2 channel - shaped pulse 90 degree (Ca on resonance)
SPNAM 3	Q3.1000	E	File name for SP3
SPOAL3	0.500	Phase alignment of freq. offset in SP3	
spoffs3 [Hz]	0	spoffs3=0	
SPW3 [W, -dBW]	21.84	-13.39	F2 channel - shaped pulse 180 degree (Ca on resonance)
SPNAM 5	Q3.1000	E	File name for SP5
SPOAL5	0.500	Phase alignment of freq. offset in SP5	
spoffs5 [Hz]	15090.93	spoffs5=bf2*(cnst21/1000000)-o2	
SPW5 [W, -dBW]	21.84	-13.39	f2 channel - shaped pulse 180 degree (C=O off resonance)
SPNAM 8	Q5tr.1000	E	File name for SP8
SPOAL8	0	Phase alignment of freq. offset in SP8	
spoffs8 [Hz]	0	spoffs8=0	
SPW8 [W, -dBW]	26.95	-14.31	F2 channel - shaped pulse 90 degree (Ca on resonance)

13C Parameter

```
(center (p14:sp3 ph1):f2 (p22 ph1):f3 )
d23
(p21 ph1):f3

(p13:sp2 ph4):f2
d0
(center (p14:sp5 ph1):f2 (p22 ph8):f3 )
d0
4u
(p14:sp3 ph1):f2
DELTA
(p14:sp5 ph1):f2
4u
(p13:sp8 ph1):f2

4u do:f1
(p26 ph7):f1
4u
p16:gp2
dl6
(p26 ph2):f1
20u cpdsl:f1 ph1
```




Channel f3

SFO3 [MHz]	50.6838141		Frequency of ch. 3
O3 [Hz, ppm]	6081.33	120.000	Frequency of ch. 3
NUC3	15N	Edit...	Nucleus for channel 3
CPDPRG 3	garp	... E	File name for cpd3
P21 [μsec]	25.000		F3 channel - 90 degree high power pulse
p22 [μsec]	50.00		F3 channel - 180 degree high power pulse
PCPD3 [μsec]	190.00		F3 channel - 90 degree pulse for decoupling sequence
PLW3 [W, dB]	135	-21.30	F3 channel - power level for pulse (default)
PLW16 [W, dB]	2.34	-3.69	F3 channel - power level for CPD/BB decoupling

Gradient channel

GPNAM 1	SMSQ10.100	... E	SMSQ10.100
GPZ1 [%]	50.00		50%
GPNAM 2	SMSQ10.100	... E	SMSQ10.100
GPZ2 [%]	40.00		40%
GPNAM 3	SMSQ10.100	... E	SMSQ10.100
GPZ3 [%]	60.00		60%
GPNAM 4	SMSQ10.100	... E	SMSQ10.100
GPZ4 [%]	30.00		30%
P16 [μsec]	1000.000		Homospoil/gradient pulse [1 msec]

15N and Gradient Parameter



Experiment

PULPROG	hncagpwg3d ... E			Current pulse program
AQ_mod	DQD			Acquisition mode
FnTYPE	traditional(planes)			nD acquisition mode for 3D etc.
FnMODE	States-TPPI States-TPPI			Acquisition mode for 2D, 3D etc.
ProjAngle [degree]	0			Angle for projection-spectroscopy
TD	1024	52	64	Size of fid
DS	32			Number of dummy scans
NS	32			Number of scans
TD0	1			Loop count for 'td0'
TDav	0			Average loop counter for nD experiments

Width

SW [ppm]	12.0161	26.0000	32.0000	Spectral width
SWH [Hz]	6009.615	1317.779	4024.467	Spectral width
IN_F [μsec]		758.85	248.48	Increment for delay
AQ [sec]	0.0851968	0.0197302	0.0079514	Acquisition time
FIDRES [Hz]	11.737530	50.683815	125.764580	Fid resolution
FW [Hz]	125000.000			Filter width

Experiment Parameter



北京大学
PEKING UNIVERSITY

上机练习HNCA



北京大学
PEKING UNIVERSITY

谢谢大家！