

Unlock Solid-State NMR: TopSolids





Jochem Struppe
Bruker BioSpin Billerica, USA
Bruker Users Meeting China, July 23th – 26th 2019



The Analytical Tool Solid-State NMR



- Solid-State NMR is often the only method of choice for solid samples when other analytical techniques cannot provide an answer:
 - Long-range order (X-ray crystallography), IR and Raman Spectroscopy and DTGA and related methods
- Uses Magic Angle Spinning (MAS) of the sample to average unwanted interactions resulting in resolved NMR signals
- Powerful method in many research areas:
 - ✓ Batteries
 - ✓ Polymers
 - ✓ Proteins
 - ✓ Pharmaceuticals



Solid-State NMR – Cliff Hanger.



- Manual probe setup (magic angle-setting, shimming)
- Thinking about, calculating & optimizing suited conditions (Hartmann-Hahn, ^1H Decoupling, etc.)
- Complicated parameter handling between different data sets
- Knowledge & expertise needed


Is Solid-state NMR - A method for experts only?

TopSolids: The Solid State NMR Expert



- The delicate setup of solid-state NMR experiments is made easy
- By reliable automated procedures
- Whatever your background and experience may be
- 2nd generation available with TopSpin 3.5.6 and TopSpin 4
- Free of charge
- Compatible with 4 mm and 3.2 mm double/ triple channel SB/ WB probes

TopSolids– Workflow




→ **TopSolids is workflow-oriented**

TopSolids – Flow Bar Modules



Different modules in flow bar design to ensure a proper workflow:

- **Preparation** (incl. advanced user interface to fetch general information in one step, import of former data possible)
- **Probe Setup** (not only for general setup, but for e.g. testing a probe after repair, etc.)
- **Applications:**
 - **Proteins** (incl. setup and optimization, up to 4D ^{13}C -detected experiments)
 - **Materials** (1D/ 2D X, HX CP, HetCor and MQMAS experiments)




TopSolids – Parameter Handling



TopSolids parameter handling is based on:

1. Parameter sets

Experiments with
pre-defined
starting values



Do not forget to run
'expinstall' for solid
systems once...

TopSolids – Parameter Handling



TopSolids parameter handling is based on:

1. Parameter sets
2. PROSOL table ('edprosol')

Enter 90° pulse
parameters for
nuclei (^1H , ^{13}C , X)

The screenshot shows the 'edprosol' window with the '90 deg. Pulses' tab selected. The window displays parameters for a 'Probe: H117468_0001 PH MASDVT 600W2 BL3.2 N/D/C/H CGR' and 'Solvent: generic'. The 'Observe' section is set to '13C' and the 'Decouple' section is set to '1H'. The 'Observe Comment' is 'Default 13C solid state observe 600 MHz' and the 'Decouple Comment' is 'Default 1H solid state decouple 600 MHz'. The '90 deg. Pulses' tab is highlighted, and a red circle around it has an arrow pointing to the text 'Enter 90° pulse parameters for nuclei (1H, 13C, X)'. Below the tabs, there is a table with columns for Nucleus, Pulse Width[μs], Power[W], and Set. The table lists parameters for 1H, 2H, 13C, 15N, and 79Br. The 13C row is highlighted.

Observe				Decouple			
Nucleus	Pulse Width[μs]	Power[W]	Set	Pulse Width[μs]	Power[W]	Set	Nucleus
1H	2.50	80.000		2.50	80.000		1H
2H	0.00	0.0000		0.00	0.0000		2H
13C	5.00	200.00		5.00	200.00		13C
15N	6.50	200.00		6.50	200.00		15N
79Br	3.80	300.00		3.80	300.00		79Br


TopSolids – Parameter Handling



TopSolids parameter handling is based on:

1. Parameter sets
2. PROSOL table ('edprosol')
3. Project-specific "parameter data bank"

- '*topsolidsPars.xml*'
- Automatically created & managed by TopSolids (not meant for manual editing)
- Stored in the project's directory
- Used to:
 - Store optimized parameters for later use
 - Import parameters from a previous project
 - Create PDF report




TopSolids – Parameter Handling



TopSolids parameter handling is based on:

1. Parameter sets
 2. PROSOL table ('edprosol')
 3. Project-specific "parameter data bank"
- Ensures reasonable values for parameters at any time – whether optimized or not yet
 - Determined parameters are stored in the parameter data bank for further usage
 - Parameters can be reviewed in form of a PDF report at any time to keep track and to assist fast publication


→ **Let's get started and create a TopSolids project...**



Reporting


No structure available.

TopSolids – Probe Setup Flow Bar




- Probe performance needs to be tested on a regular basis to ensure best results
- Though routine work, it can be time consuming - even for an NMR expert
- The ***Probe Setup*** module offers a high level of convenience to:
 - ✓ Setting the magic angle
 - ✓ Optimizing ^1H and ^{13}C hard pulses on standard samples
 - ✓ Shimming the probe and referencing the magnetic field
 - ✓ With least user interaction to overcome the necessity for extensive prior knowledge
 - ✓ Only two samples needed: KBr, adamantane


Probe Setup – Set Magic Angle



- On KBr (^{79}Br) at 5 kHz MAS TopSolids:
 - ✓ Checks for fixed frequency H/C/N probes to adjust respective parameter sets automatically
 - ✓ Detects Cover Module 2, which can be used for MA adjustment at the spectrometer directly – if wanted
 - ✓ Records a reference spectrum and evaluates the quality of the MA setting
 - ✓ Performs an angle-stability test for SB probes
 - ✓ Checks final quality



Probe Setup – Set Magic Angle using the *iProbe*



- Using an iProbe, TopSolids:
 - ✓ Automatically tunes & matches the probe
 - ✓ Automatically adjusts the angle setting
 - ✓ Records a reference spectrum and evaluates the quality of the MA setting
 - No manual interactions required



Probe Setup – Set Magic Angle using the *iProbe*



Acquire Process **>click<** Applications Manage

Exit Probe Setup **Set Magic Angle** 1H Hard Pulse 13C Hard Pulse Shimming & Referencing Go On

MAS display

Dataset
Select a dataset for Magic Angle adjustment
Create a new dataset for Magic Angle adjustment

Start adjustment
Start automatic Magic Angle adjustment
Start manual Magic Angle adjustment

Stop adjustment
Stop Magic Angle adjustment

Adjust
Magic Angle
Reset

Debug
Minimum position: 141500 Maximum position: 147500 Current position: 147035 Default Magic Angle position: 144635

Acquisition window

title information updated frequently

shows progress of the automated angle setting

shows active commands

SOL_ANGLE_SETTI

Please insert and rotate a KBr sample at 5 kHz MAS

As soon as the MAS rate is stable, click 'OK' to continue.

'Cancel' exits the setup process.

OK Cancel

topsolids 15 1 C:\NMRDATA\data\MAS40\nmr

Standard Setup
Magic Angle setting
acquisition time: 40.96 ms
number of time domain points: 8192
LW40.41 (line width at half height center peak): 116.18 Hz
LW40.42 (line width at half height 5th left spinning side band): 148.39 Hz

Acquisition
PULPROG = C.90
NUC1 = 79Br
SU = 996.99
SWM = 100000
TD = 8192
SCANS = 16/16

0.01 0.02 0.03

topsolids 15 1 C:\NMRDATA\data\MAS40\nmr SOL_ANGLE_SETTI in progress

Network C:\ML>assist Login

MAS display

Main Magic Angle Monitoring Recording Configuration Log Help

Dataset

Select a dataset for Magic Angle adjustment

Create a new dataset for Magic Angle adjustment

Start adjustment

Start automatic Magic Angle adjustment

Start manual Magic Angle adjustment

Stop adjustment

Stop Magic Angle adjustment

Adjust

Magic Angle

<<<< <<< << << < Reset > >> >>> >>>>

Debug

Minimum position: 141500 Maximum position: 147500 Current position: 147035 Default Magic Angle position: 144635

MAS Unit State: ☒ Rotation Running | Spinrate: 5001 Hz | System Pressure: 6780 mbar | Probe Temperature: Initializing | Recording

SPECTRUM PROCPARS ACQPARS TITLE PULSEPROG PEAKS INTEGRALS SAMPLE STRUCTURE PLOT FID **ACQU**

topsolids 15 1 C:\NMRDATA\data\MAS40\nmr

Standard Setup
Magic Angle setting
acquisition time: 40.96 ms
number of time domain points: 8192
LYWH-1 (line width at half height center peak): 116.18 Hz
LYWH-2 (line width at half height 5th left spinning side band): 148.39 Hz

Acquisition

PULPROG = C.90
NUC1 = 79Br
SV = 996.99
SWH = 100000
TD = 8192
SCANS = 16/16


SOL_ANGLE_SETT

Please insert and rotate a KBr sample at 5 kHz MAS.
As soon as the MAS rate is stable, click 'OK' to continue.
'Cancel' exits the setup process.


0.01 0.02 0.03 [s]

topsolids 15 1 C:\NMRDATA\data\MAS40\nmr SOL_ANGLE_SETT: in progress


Probe Setup – Hard Pulse Optimization




- Determining the 90 degree pulses for ^1H & ^{13}C is:
 - ✓ Done automatically
 - ✓ On adamantane at 10 kHz MAS
 - ✓ Using the script routine '*pulsecal*' (known from solution NMR, adjusted for TopSolids)
 - ✓ Using PROSOL values as starting parameters or data bank entries IF present (= imported)
 - ✓ Taking ~2 minutes per nucleus




Probe Setup – Shimming & Referencing




- Shimming essential to ensure for well-resolved lines
- Done on adamantane at 10 kHz MAS
- Automated Z shimming routine available (2-5 min depending on starting values)
- Shim quality determined automatically & compared to starting values
- B_0 field is automatically referenced to neat TMS



Note: Z shim values will be updated by the routine






TopSolids – Applications Module



- Can be opened as soon as 'Shimming & Referencing' in the 'Probe Setup' module has been performed ✓
- Dedicated workflows for:
 - ✓ Protein samples
 - ✓ Non-protein samples ('Materials')



TopSolids – Protein Applications



Whole setup divided into three modules:

1. 'Standard Sample Setup'


- Should be performed if the sample of interest is giving poor signal to noise

2. 'Sample of Interest Setup'


- (Re-) optimize parameters on your sample of interest
- Optional optimizations: spectral widths, carrier frequencies, ^1H decoupling during acquisition

3. 'nD Experiments'



- Actual experiments of interest
- 1-click-optimization possible → does no longer rely on prior optimization




Protein Applications – Standard Sample Setup



- Standard sample: [^{13}C , ^{15}N]-glycine
- For NCO transfer: any [^{13}C , ^{15}N]-peptide (e.g. fMLF) or standard protein (e.g. SH3, GB1, ubiquitin)
- Cross polarization (CP) transfers:
 - 'Parameter Check' based on optimized hard pulses & current MAS rate
 - For MAS III units the rate is detected automatically
 - Else: put in manually





Protein Applications – Standard Sample Setup




Based on 'Parameter Check'
TopSolids suggests a suitable
Hartmann-Hahn condition:


- ✓ Checks power limits
- ✓ Checks for unwanted conditions
- ✓ Free choice of shape
- ✓ Takes into account ^1H decoupling during double CP






Protein Applications – Sample of Interest Setup



- Comprises optimizations of all parameters needed for multidimensional experiments available in TopSolids
- General setup of optimizations equals the Standard Sample Setup
- If the 'Standard Sample Setup' module has not been executed before, starting values are taken from PROSOL (+ parameter data bank, if parameters are available)



Protein Applications – Sample of Interest Setup



- The smaller a spectral width in an indirect dimension, the shorter the experimental measurement time of an nD experiment
- TopSolids optimizes spectral widths & carrier frequencies fully automated - based on CP experiments

Protein Applications – nD Experiments Module



Navigation bar: Acquire | Process | Analyse | Applications | Manage

Exit | Proteins - nD Experiments (13C) | Intraresidue Backbone Assignment | Interresidue Backbone Assignment | Distance Measurements | Go On

2D Experiments ▶	PDSD/ DARR (short C-C mixing)
3D Experiments ▶	NCa
	NCaCx

2D Experiments ▶	NCaCx
3D Experiments ▶	

2D Experiments ▶	PDSD/ DARR (long C-C mixing)
3D Experiments ▶	NCO
4D Experiments ▶	NCOCx

2D Experiments ▶	NCOCx
3D Experiments ▶	CaNCOCx
4D Experiments ▶	CONCaCx

2D Experiments ▶	CaNCOCx
3D Experiments ▶	CONCaCx
4D Experiments ▶	

PDSD/ DARR (very long C-C mixing)
ChhC
NhhC

- 2D to 4D experiments for protein backbone assignment and distance measurements

nD Experiments Module – '1-Click-Optimization'



- Do I need to go through all these optimizations if I only want to record one specific nD experiment?
- Which parameters do I need for recording a 3D NCOCx experiment?
- ✓ Skip the Sample Setup modules & start directly with the 'nD Experiments' using '1-click-optimization'!
- ✓ Select an nD experiment of interest in the 'nD Experiments' module (e.g. 3D NCOCx)
- ✓ TopSolids automatically detects missing parameters & initiates the respective optimization directly using least user interaction
- ✓ Once all parameters have been optimized, it will go back to the nD experiment to support its set up

nD Experiments Module – '1-Click-Optimization'



Spectrum ProcPars AcqPars Title PulseProg Peaks Integrals Sample Structure Plot Fid Acqu

Backbone Assignment

3D NCOCx correlation experiment (CC transfer via PDSD)
on biological sample of interest

CO-Cx mixing time: 30 ms

MAS rate: 10.0 kHz

number of scans: 48

relaxation rate: 2.0 s

acquisition time F3(Cx): 22.53 ms

acquisition time F2 (CO): 5.70 ms

acquisition time F1 (N): 5.60 ms

number of time domain points F3 (Cx): 2048


number of time domain points F2 (CO): 38

number of time domain points F1 (N): 28

total experimental time: 1 d 5 h 46 min 42 sec

experiment in queue

experiment finished





- 3D NCOCx recorded with 1-click-optimized parameters
- Title includes information about acquisition parameters
- Processing was done fully automated by TopSolids after acquisition was finished
- TopSolids automatically queues experiments in the spooler when acquisition is active

Processing


TopSolids – Materials Applications




- Any X nucleus of choice the mounted probe can handle
 - TopSolids automatically detects, if the chosen X nucleus is:
 - A spin $1/2$ nucleus
 - A quadrupolar nucleus
 - Compatible with the chosen probe
- Optimizations available for:
 - ✓ X 90 degree pulse
 - ✓ HX Cross Polarization (CP)
 - ✓ 2D HetCor
 - ✓ 1D/ 2D MQMAS






Materials – Direct Excitation Module




Materials – Direct Excitation Module



- Guides through basic setup routine of 90 degree pulse calibration on a standard samples and/ or the sample of interest
- X 90 degree pulse optimization either as a high or a low power pulse, depending on the X nucleus of choice
- (High/ low power) ^1H decoupling optional




Materials – CP-based Experiments




- Based on successful performance of X 90 degree pulse calibration
 - TopSolids automatically starts the '1-click-optimization' if missing
- CP transfer can be optimized on a standard samples and/ or the sample of interest
- CP conditions differ for spin $\frac{1}{2}$ and quadrupolar nuclei
 - TopSolids is taking care to guarantee for best transfer efficiency

Materials – CP-based Experiments




- Based on the H-X CP parameters a 2D H-X HetCor spectrum can be set up and recorded
- Spectrum automatically processed and phased




Processing


Materials – MQMAS Experiments



- Stand-alone module in the Materials Applications
- Guides stepwise through the setup to run a 2D MQMAS experiment:
 1. X 90 soft pulse (if not known already) on standard sample and/ or sample of interest
 2. Specific MQMAS excitation & reconversion pulse lengths for optimum sensitivity
 3. Minimum spectral width to reduce 2D measurement when used as indirect dimension




Materials – MQMAS Experiments



4. Stepwise 2D setup where TopSolids:


- Takes care that a multiple of the phase cycle is used as number of scans
- Automatically synchronizes the optimized spectral width with the MAS rate
- Applies basic processing and shearing
- User interaction only needed for entering number of scans & relaxation delay




TopSolids – Navigation



- Whenever TopSolids shall be restarted, open an ExpNo of the respective project directory and start TopSolids
- TopSolids recalls the module that has been used latest
- Navigation between the different Application modules is done easily...





TopSolids – Features at a Glance



- Intelligent handling of project and parameters (incl. import/ report)
- Probe setup: within less than 30 minutes possible
 - Automated quality evaluation of magic angle adjustment and shimming
- Protein Applications: state-of-the-art pulse program library up to 4D experiments
 - “1-click-optimization” – TopSolids detects missing parameters to initialize their optimization
 - Automated optional optimizations of spectral widths, carrier frequencies and decoupling parameters
- Materials Applications: suited for any X nucleus the hardware can measure, incl. quadrupolar nuclei
 - High/ low power X 90 degree pulses
 - HX cross polarization (CP) transfer & 2D HetCor
 - 1D/ 2D MQMAS experiments
- Queuing of multidimensional experiments
- Automated processing up to 3D, including phase corrections (& shearing for MQMAS experiments)

TopSolids – Benefits at a Glance



- Quality output
- Minimized risk of damage
- Efficient use of spectrometer time
- Uniform workflows
- Keep track via fast reporting

... whatever your background and experience will be

→ TopSolids: Your key to unlock solid-state NMR !

Further Information



- Do you want to get in contact with us?
 - solids@bruker.com
- Do you want to get further details?
 - Product page at www.bruker.com

The screenshot shows the Bruker website's navigation bar with links for Products, Applications, Service, News, Events, About us, and Login. A search bar is on the right. The main banner for TopSolids features a blue background with the text "TopSolids" and a "Contact" button. Below the banner is a breadcrumb trail: Home - Products - Magnetic Resonance - NMR - Software - TopSolids - Overview. A secondary navigation bar includes tabs for Overview, Technical Details, Application, and Learn More, along with a contact icon. The main content area is titled "TopSolids™ - Accelerated Workflow in Solid-State NMR" and contains three columns of text describing the software's capabilities.

TopSolids™ - Accelerated Workflow in Solid-State NMR


Automated spectrometer adjustment together with a comprehensive range of pre-defined protocols guarantees optimum spectral quality and high productivity for both expert and novice users

TopSolids™ is an intuitive, workflow-based software package for solid-state NMR that is especially designed to give an easy access into the delicate setup of probe hardware and difficult experiments dedicated to insoluble protein systems and solid samples, e.g. from pharmacy, catalysis, batteries or material

TopSolids™ ensures robust, secure spectrometer and probe setup and guarantees unsurpassed spectral quality by combining optimized acquisition parameters with fine tuning on the actual research sample.

With TopSolids™ Bruker offers a high level of convenience to use one of the most important

谢谢



Innovation with Integrity

结束